

Intestinal epithelial cells – primary cell culture

1. Flush the contents of the intestine with Ca^{2+} - and Mg^{2+} -free Hanks' balanced salt solution (HBSS) containing 2% glucose, 25 ng of amphotericin B per ml, 100 U of penicillin per ml, and 100 μg of streptomycin per ml.
2. Splice the intestine into small pieces and incubate for 15 min at 22°C on a shaker platform in Ca^{2+} - and Mg^{2+} -free HBSS containing 5mM EDTA, 2% bovine serum albumin, and 0.2 mg of soybean trypsin inhibitor per ml.
3. Take out the supernatant and wash with (DMEM), 100 U of penicillin per ml, 100 μg of streptomycin per ml, and 5% fetal bovine serum (FBS)
4. Cells were cultured in 24-well plates, at a seeding density of approximately 2×10^6 cells/well
5. One hour before plating cells, culture surfaces were coated with 40 μl of Matrigel (BD Biosciences) per cm^2 diluted 1:2 in phenol-red-free DMEM (Sigma).
6. Epithelial cells were cultured in epithelial cell medium (ECM) containing equal volumes of phenol-red-free DMEM and Ham's F-12 medium (Biowhittaker) with the following additives:
 - 5 μg of insulin (Sigma) per ml,
 - 5×10^{-8} M dexamethasone (Sigma),
 - 60 nM selenium (Sigma),
 - 5 μg of transferrin (Sigma) per ml,
 - 5×10^{-8} M triiodothyronine (Sigma),
 - 10 ng of epidermal growth factor (Sigma) per ml,
 - 20 mM HEPES,
 - 2 mM glutamine,
 - 100 U of penicillin per ml,
 - 100 μg of streptomycin per ml,
 - 0.2% D-glucose,
 - 2% FBS.
7. Cells were cultured in 5% CO_2 at 37°C with periodic supplementation of medium to maintain a volume of 2 ml per well.