

# Immunoprecipitation Protocol

## Reagents Needed:

Immunoprecipitation (IP) lysis buffer  
Protease Inhibitors (Calbiochem Cat#539131)  
Primary Antibodies made in Rabbit  
Normal IgG, negative control (Rabbit IgG- Bethyl Cat. No. P120-101)  
Protein A Sepharose Beads (Amersham Cat# 17-0780-01)  
Cell Lysate  
Sample Buffer

## IP lysis Buffer

12.5 ml 1M NaCl (250mM)  
2.5 ml 1M Tris (50mM)  
500ul 0.5M EDTA (5mM)  
2.5ml 10% NP-40  
32 ml dH<sub>2</sub>O

## Protein A Beads

Resuspend 400 mg of Protein A beads in 10 ml of distilled H<sub>2</sub>O. Mix well to resuspend. Spin at 250 rpm for 5 minutes. Wash 3X in 10 ml IP Lysis buffer. Resuspend to 10 ml with IP lysis buffer for a 20% solution. Use 100 µl per IP reaction.

## 4X Sample Buffer (Store at 4 C)

Glycerol - 4.0 g  
Tris Base - 0.68 g  
Tris HCL - 0.67 g  
LDS - 0.80 g  
EDTA - 6 mg  
Brilliant Blue G250 - 2.5 mg  
Phenol red - 2.5 mg

## 1X Sample Buffer

4X sample buffer - 150 µl  
1M DTT - 60 µl  
Distilled water - 390 µl  
Make **fresh** for each use.

## Procedure:

1. Place 500 µl of the pooled cell lysate (1-3 mg/ml) into a 1.5 ml micro-centrifuge tube.
2. To this tube add 2 to 10 µg of the primary antibody (If using neat sera or an IgG fraction such as Protein-A purified antibody, larger amounts are likely to be required. For best results, optimal amounts of antibody should be empirically defined.)
3. To a negative control reaction, add an equivalent amount of normal rabbit IgG.

4. Add 100  $\mu$ l of a 20% Protein A suspension.(Amersham Biosciences, Cat# 17-0780-01) to the mixture of antibody and cell lysate. Rotate the immunoprecipitation reactions (end-to-end) for 3 hours at room temperature or overnight at 4 C.
5. Centrifuge (200 x g; 5 minutes) to pellet the complex.
6. Remove the supernatant and add 500  $\mu$ l cold cell lysis buffer. Centrifuge (200 x g; 5 minutes).
7. Repeat wash step 6 twice more. After each centrifugation remove as much of the supernatant as possible.
8. After removing the supernatant from the third wash, add 40  $\mu$ l of freshly prepared 1X sample buffer to each tube and heat at 90 C for 5 minutes.
9. Continue with electrophoresis and immunoblotting as described under western blotting protocol. Load 8 to 16  $\mu$ l (20 to 40% of the IP reaction) to a polyacrylamide gel.

**Note:** For optimal results, complete reduction of the sample is required. We recommend the use of 0.1 M DTT in SDS-PAGE sample buffer and immediately heating samples, loading and running gels